


Product Name	Recombinant Human CD47	SDS-PAGE gel	
Synonym(s)	IAP; MER6; OA3	1	2
Quantity	100 µg, 1 mg	kDa	M
Species	Human	116	
Class Type	Recombinant Protein	66.2	
Molecular weight	48-60 kDa (under reducing conditions), 41 kDa (predicted)	45.0	
Purity	>95% by SEC-HPLC, >90% by SDS-PAGE	35.0	
Tag	C-terminal human IgG1 Fc tag	25.0	
Expression Source	HEK 293 cells	18.4	
GenBank Accession #	NP_942088.1 (Met1-Pro139)	14.4	
Application	Western Blotting, Functional Studies		
Formulation	Lyophilized from sterile PBS, pH 7.4.		
Storage and Stability	Stable for 12 months at -80°C, Avoid freeze/thaw cycles	1 – MW Marker 2 – CD47	
Description	CD47 contains 1 Ig-like V-type (immunoglobulin-like) domain and is a receptor for the C-terminal cell binding domain of thrombospondin. It may play a role in membrane transport and signal transduction. CD47 is also a membrane protein, which is involved in the increase in intracellular calcium concentration that occurs upon cell adhesion to extracellular matrix. It is very broadly distributed on normal adult tissues, as well as ovarian tumors, being especially abundant in some epithelia and the brain. CD47 may play a role in membrane transport and/or integrin dependent signal transduction. It may prevent premature elimination of red blood cells. It also may be involved in membrane permeability changes induced following virus infection. By acting as an adhesion receptor for THBS1 on platelets, CD47 plays a role in both cell adhesion and in the modulation of integrins. It also plays an important role in memory formation and synaptic plasticity in the hippocampus.		
Reference	Brown EJ, et al. (2001) Integrin-associated protein (CD47) and its ligands. Trends Cell Biol. 11(3): 130-5. Oldenborg PA. (2004) Role of CD47 in erythroid cells and in autoimmunity. Leuk Lymphoma. 45(7): 1319-27. Kaczorowski DJ, et al. (2007) Targeting CD47: NO limit on therapeutic potential. Circ Res. 100(5): 602-3.		

Related products:

Product Name	Catalog #	Size
Recombinant Human PD-L1	237351	100 ug
Recombinant Human CD40	232340	100 ug
Recombinant Human CD40L	2323405	100 ug
Recombinant Human SIRP-alpha	2374772	100 ug

Recombinant Human CD155	2323155	100 ug
Recombinant Human GITRL	2344875	100 ug
Recombinant Human OX40	236940	100 ug
Recombinant Human OX40L	2369405	100 ug
DNA Polymerase Theta Activity Assay Kit	362101	96 reactions
DNA Polymerase Theta-N-Helicase Domain	7657643	20 ug, 100 ug
DNA Polymerase Theta-C terminal Domain	7657283	20 ug, 100 ug, 1 mg
DNA Polymerase Theta Full Length protein	7657385	10 ug, 50 ug
T7 RNA polymerase	777627	5000 U, 25000U, 100000U
T7 High Yield RNA Synthesis Kit	777627-RK	25 rxns, 50 rxns, 100 rxns
Recombinant Mouse Leukemia Inhibitory Factor (mLIF)	11-0001	10 ug, 100 ug
Recombinant Human LIF	12-0002	10 ug, 100 ug, 1 mg
Recombinant Human LIF, Animal-Free	12-0002AFR	10 ug, 100 ug, 1 mg
Recombinant Human FGF-basic, Carrier-free	12-0005CFR	50 ug, 100 ug, 500 ug, 1 mg
Kras Wild Type (WT), GST-tag	5727-4121G	50 µg, 100 µg
Kras WT, GST-tag, GDP Loaded	5727-WTG-G	50 µg, 100 µg
Kras WT, GST-tag, GppNHp loaded	5727-WTG-GP	50 µg, 100 µg
Kras G12C, GST-tag	5727-4122G	50 µg, 100 µg
Kras G12C, GST-tag, GDP Loaded	5727-4122G -G	50 µg, 100 µg
Kras G12C, GST-tag, GppNHp loaded	5727-4122G -GP	50 µg, 100 µg
Kras G12C Nucleotide Exchange Assay Kit	5727-4122NK	384 reactions
Kras G12C – cRAF Binding Assay Kit	5727-4122BK	384 reactions
Human RBD-RAF1, N-His tag, C-FLAG tag	7237231	50 µg, 100 µg
Human SOS1, No Tag	7671	50 µg, 100 µg
Human SOS1, His-Avi-Tag	7671HA	50 µg, 100 µg

This product is for research use only and not for diagnostic or therapeutic use.

Aurora Biolabs LLC, San Diego, CA 92121, USA; www.aurorabiolabs.com;

order@aurorabiolabs.com; info@aurorabiolabs.com Tel: 858-215-4510 or 858-453-5700 Fax: 855-898-3979